

MyBabyMonitor



Team Members:


Benedetta Boccardi – 292175

Martina Benvenuto – 291327


Nicola Battaglini – 287468

Domenico Ficili – 279931

The idea



Causes




Solution

A continuous monitoring system for evaluation of vital signal monitoring can help *reduce* premature infant mortality rate.



Features



Software characteristics



Integration

External services can be integrated through APIs



Modularity and Scalability


- Registration of multiple devices and users
 - Simultaneously handling of data coming from different devices
- Simple integration of new functionalities



Microservices Architecture

Services run independently


System Overview



Use Case Diagram: proposed

Legend:


- — REST Web Services (provider)
- ⌋ — REST Web Services (consumer)
- - - MQTT Communication
- (1) Publisher
- (2) Subscriber



Use Case Diagram: modified

Legend:

- — REST Web Services (provider)
- ⌋ — REST Web Services (consumer)
- - - MQTT Communication
- (1) Publisher
- (2) Subscriber



Device



Sensors




Heartrate sensor



Pulsossimeter: oxygen level



Microphone



MIMU sensor: seizure



simulated in .txt files

Device

1

MQTT: Publish

- Publish messages relative to the measurements using SenML format to the microservices with the topic "MyBabyMonitor/**Service/DeviceID** "
- Publish the event with the topics "MyBabyMonitor/Alarm/**DeviceID**" and "MyBabyMonitor/Alarm/**DeviceID/Duration**"

2

MQTT: Subscribe

Subscribes to the topic "MyBabyMonitor/Event/**DeviceID**" to receive events from the microservices

3

REST


The Blood Level Oxygen and Heart Rate data are sent to the Thingspeak personal channel using REST protocol

SenML format



```
apneaMessage={
  "bn": deviceID,
  "e": [
    {
      "n": "respiration",
      "u": "",
      "t": time.time(),
      "v": resp
    },
    {
      "n": "oxylevel",
      "u": "",
      "t": time.time(),
      "v": oxy
    }
  ]
}
```

Catalogs




Main Catalog

- Information of registered users and devices
- System configuration data
- List of online available services



Seizure Catalog

- Seizure events for each user



Apnea Catalog

- Apnea events for each user



Monitoring Catalog

- Apnea and Seizure events data for each user

Catalog Manager



```
{
  "projectOwner": "Gruppo 6",
  "projectName": "MyBabyMonitor",
  "lastUpdate": "2022-07-01-10:43:32",
  "lastChildID": 14,
  "thingSpeak": {
    "url": "https://api.thingspeak.com/update.json",
    "userAPIKey": "PNZCK58R84TT8UKE"
  },
  "telegram": { ...
},
  "urlCatalog": "https://host.docker.internal:8080/",
  "mosquitto": { ...
},
  "onlineServices": [],
  "childrenList": [ ...
],
  "devicesList": [
    {
      "deviceID": 6,
      "childID": 5,
      "sensorsList": [
        {
          "sensorID": "6-1147",
          "sensorName": "Microphone",
          "measureType": [
            "Sound"
          ],
          "availableServices": [
            "MQTT"
          ],
          "servicesDetails": [
            {
              "serviceType": "MQTT",
              "serviceIP": "mqtt.eclipse.org",
              "topic": [
                "MyBabyMonitor/apnea/sound"
              ]
            }
          ]
        }
      ],
      "lastUpdate": ""
    }
  ]
}
```

Central unit of the system

Main functions:

- Setting services online
- Obtaining status info about services
- Updating and removing services
- Working as interface between the database and the rest of the System.
- Managing users' registration and account settings


Using its functions, it can be retrieved:

- Microservices URLs
- Thingspeak configuration
- Broker configuration (Mosquitto)
- MQTT topics used by other services


The main catalog, which serve as database, is exposed to microservices through REST

Apnea

Apnea microservice detects an event when for 5 seconds the following conditions occur:



Using **REST** interface, the service communicates to the server that is online and gets information to configure Mosquitto.




MQTT protocol is used:


- to receive messages about sensor measurements from the Raspberry device
- sort the data according to the origin using DeviceID
- to communicate to the Raspberry if an event has been detected

Seizure

Seizure microservice detects an event when for 5 seconds the following conditions occur:



Using **REST** interface, the service communicates to the server that is online and gets information to configure Mosquitto.



MQTT protocol is used:

- to receive the messages about sensor measurements from the Raspberry device
- sort the data according to the origin using DeviceID
- to communicate to the Raspberry if an event has been detected

Monitoring

Monitoring microservice handles visualization of:


- real-time acquired data
- stored data events

Using **REST** interface, the service communicates to the server that is online and gets information to configure Mosquitto.


MQTT protocol is used to receive alarm messages that are stored in the MonitoringDB.

It subscribes to the topic:
"MyBabyMonitor/Alarm/#"


The information stored in the MonitoringDB can be visualized through the website




Monitoring: Website



Home



Website is built using Flask and Bootstrap Framework




Apnea/Seizure




Website retrieves data of events and shows them in two tables

Monitoring: Thingspeak




Thingspeak collects data received from the user's device, store them and creates graphs in the personal channel



The Monitoring service retrieves graphs from the user's Thingspeak channel and shows them in the personal webpage


Telegram Bot: User Interface

MyBabyMonitorBot serves as User Interface to register and manage users' account




Telegram Bot: Notification

MQTT protocol is used to receive alarm messages and send a notification to the user using Telegram Bot.



It subscribes to the topic:
"MyBabyMonitor/Alarm/#"


Notification alarm



Docker



Containers



Docker containers are used to:

- Deploy system microservices
- Run each service in an independent container
- Distribute the workload and speed up the system implementation

Conclusion

STRENGTHS

- Real-time monitoring of baby's vital parameters with embedded graphics
- Detection algorithm and notification system of apnea and seizure events
- Microservices architecture
- Modular, flexible and scalable system

S

WEAKNESSES

- Data sent and stored in ciphertext on the database
 - System deployed to run on local machine

W

- Implementing cryptographics and pseudonimization algorithms to improve security of the system
- Deploying the system to run on online server

O

OPPORTUNITIES

- System vulnerability to malicious attacks

T

THREATS

Bibliography

- [1] S. Shamsir, O. Hassan and S. K. Islam, "Smart Infant-Monitoring System with Machine Learning Model to Detect Physiological Activities and Ambient Conditions," 2020 IEEE International Instrumentation and Measurement Technology Conference (I2MTC), 2020, pp. 1-6, doi: 10.1109/I2MTC43012.2020.9129295.
- [2] Gravem, D., Singh, M., Chen, C., Rich, J., Vaughan, J., Goldberg, K., Waffarn, F., Chou, P., Cooper, D., Reinkensmeyer, D., and Patterson, D. (May 14, 2012). "Assessment of Infant Movement With a Compact Wireless Accelerometer System." ASME. J. Med. Devices. June 2012; 6(2): 021013. <https://doi.org/10.1115/1.4006129>
- [3] Breathing and heart rates in unwell children site: <https://www.northlondonpaediatrician.co.uk/wp-content/uploads/2014/07/NLP-heart-and-resp-rates.pdf>

**Thanks for the
attention!**



Any question?